


Nanodyne Replacement Illuminator for B & L StereoZoom Illuminator: Included Items

PN 10238 B&L Illuminator Assy



Illuminator installed on microscope - microscope not included.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	3/14/2019
PN 10518 B & L Illuminator System	REV 3
SHEET 1	OF 1