

Nanodyne Replacement Illuminator for Diagnostic Instruments TLB D4.1 Microscope Base: Included Items

Additional Items Included But Not Shown:

- PN 10709 Hex Key 3/32 Short Arm (for 4-40 screws)
- PN 11564 Hex Key 7/64 Long Arm (for 6-32 screws)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.


Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.



Han-Seung Yang		1/21/2020
PN 11559 TLB D4.1 illuminator System	REV	1
SHEET 1		OF 1