


Nanodyne Replacement Illuminator for Fisher Scientific Inverted Microscope - Included Items

PN 10334 Olympus BX40 Series Illuminator Assy




The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

PN 11521 P160KNP Pot-Cable Assy - shortened 5mm with Tab (mod with 34 inch long cable)



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		10/25/2019
PN 11293 Fisher Scientific Illuminator System		REV 3
SHEET 1		OF 1