


Nanodyne Replacement Illuminator for Labomed CXR 3 Microscope: Included Items

Tools not included but required:
 PN 10337 Hex Key 2mm Short Arm for M3 BHCS
 Phillips Screwdriver
 11 mm socket
 Misc tools for removing panels and original illuminator from microscope.



3M VHB Tape to secure wires.



PN 10733 Power Supply & PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		12/20/2023
PN 11831 Labomed CXR 3 Illuminator System		REV 2
SHEET 1		OF 1