


Nanodyne Replacement Illuminator for Leica DMLS Microscope: Included Items

Additional Items Included But Not Shown:

- PN 10456 1.5mm Hex Key (for pot adjustment knob)
- PN 10337 2mm Hex Key (for M3 mounting screws)
- PN 10444 Tool to remove OEM lamp assembly



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		11/30/2020
PN 10499 Leica DMLS Full Illuminator System		REV 4
SHEET 1		OF 1