

# Nanodyne Replacement Illuminator for Leica DMLB 100W Microscope - Included Items

## Additional Items Included But Not Shown:

- PN 10456 Hex Key 1.5mm (for pot knob)
- PN 10747 3mm Hex Key (for set screw holding lamp housing)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.



PN 10736 Rubber plug to block unused AC power receptacle.




Tape to secure wires.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.


Wayne Bonin 11/25/2019

PN 10746 Leica DMLB 100W Illuminator System REV 6

SHEET 1 OF 1



PN 10589 Leica DMLB 100W Illuminator Assy



PN 10497 Pot-Cable Assy (16 inch long)

Pot Knob from DMLB 100W (save from original microscope)