


Nanodyne Replacement Illuminator for Leica DM1000 Microscope: Included Items

Additional Items Included But Not Shown:

- PN 10457 Hex Key 2.5mm Short Arm for M3 SHCS
- PN 10337 Hex Key 2mm Short Arm for M3 BHCS



PN 10736 Rubber plug to block unused AC power receptacle.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



Tape to secure wires.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang 11/23/2020	
PN 11819 Leica DM1000 Illuminator System	REV 1
SHEET 1 OF 1	