

Nanodyne Replacement Illuminator for Leica DMIL Microscope Illuminator (35W OEM Lamp) - High Power: Included Items

Additional Items Included But Not Shown:

PN 12191 Leica DMIL  
Illuminator Assy  
(35W OEM lamp)  
High Power

PN 10709 Hex Key 3/32 inch for set screw

PN 10761  
SHCS  
4-40 x 3/16



PN 11163 Power Supply - 5V 2.1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10736 Rubber plug to block unused AC power receptacle.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	3/29/2022
PN 12192 Leica DMIL Illuminator System (35W OEM lamp) - High Power	REV 1