

Nanodyne Replacement Illuminator for Leica DMIRE2 / DMIRB Microscope: Included Items

Also Included but not shown:
 PN 10337 Hex Key 2mm Short Arm (for M3 BHCS)
 PN 10457 Hex Key 2.5mm Short Arm (for M3 SHCS)

PN 12212 Leica DMIRE2 Illuminator Assy

PN 11531 Power Cable Jack to Wires
 (6 in cable)

Leica DMIRE2 Microscope OEM bracket

Leica DMIRE2 Microscope frame

PN 11170 Pot-Cable Assy
 (detabbed, full length shaft)
 45 in cable



3M VHB Tape to secure wires.



PN 12434 Power Supply US 15V36W
 Barrel Plug 2.5mm x 5.5mm x 9.5mm
 and PN 10734 Power Cable 1.35mm
 ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	10/30/2023
PN 12213 Leica DMIRE2 Illuminator System	REV 1
SHEET 1 OF 1	