

Nanodyne Replacement Illuminator for Leitz Wetzlar SM-LUX-POL Microscope - Included Items

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch (for 8-32 set screw)

The illuminator mounted adjustment potentiometer is standard on this unit. Ask about a remote adjustment pot if desired.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin		12/9/2019	
PN 11218 Leitz Wetzlar SM-LUX-POL illuminator System			REV 2
SHEET 1		OF	1