

Nanodyne Replacement Illuminator for Leitz Wetzlar SM-LUX-POL Microscope - Included Items

PN 10983 T2PI
Main Assy R3
Side Adjust with Lens

PN 10467 Set Screw
8-32 Brass Tip

PN 11202 Leitz SM-LUX
Bottom illuminator Adapter

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch (for 8-32 set screw)

The illuminator mounted adjustment potentiometer is standard on this unit. Ask about a remote adjustment pot if desired.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	12/9/2019
PN 11218 Leitz Wetzlar SM-LUX-POL illuminator System	REV 2
SHEET 1 OF 1	