

Nanodyne Replacement Illuminator (High Power) for Leitz Wetzlar SM LUX POL Microscope - Included Items

PN 11145 T2PIHP3  
Illuminator Assy with  
Side Adjust

PN 10467 Set Screw  
8-32 Brass Tip

PN 11202 Leitz Wetzlar  
SM-LUX POL Adapter

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch (for securing illuminator assembly to adapter)

\*Note that the illuminator mounted adjustment potentiometer is standard on this unit. Ask about a remote adjustment pot if desired.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device. (USB ports may not provide full output power for high power illuminators)

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

|   |  |           |   |
|---|--|-----------|---|
| Wayne Bonin   |  | 12/9/2019 |   |
| PN 11220 Leitz Wetzlar SM LUX POL 100W illuminator System |  | REV       | 2 |
| SHEET 1   |  | OF        | 1 |