

**Nanodyne Replacement Illuminator for Leitz Laborlux Microscope (Bottom) - Included Items**

**Additional Items Included But Not Shown:**

PN 10490 Hex Key 5/64 inch for set screw  
 PN 10747 Hex Key 3mm for M4 screws

PN 11541 Leitz Laborlux illuminator  
 Assy (Bottom) - Side Adjust

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.




Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.

Han-Seung Yang 10/13/2020	
PN 11542 Leitz Laborlux illuminator System (Bottom) - Side Adjust	REV 1
SHEET 1 OF 1	

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10467 Set Screw  
 8-32 Brass Tip

PN 11500 Leitz Laborlux  
 Adapter Plate (Bottom)

PN 11504 SHCS M4 x 0.7 x 20mm  
 (2 pcs)