


**Nanodyne Replacement Illuminator for Leitz Ortholux POL Microscope Illuminator: Included Items**



**Additional Items Included But Not Shown:**

PN 10490 Hex Key 5/64 inch for set screw



PN 11163 Power Supply 5V 2.1A and  
PN 10734 Power Cable 1.35mm  
ID/3.5mm OD x USB A 6 ft.



PN 11795 Remote Control Box Assy  
**Available on request**

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	1/5/2021
PN 11908 Leitz Ortholux POL Illuminator System	REV 1
SHEET 1 OF 1	