


**Nanodyne Replacement Illuminator for Leitz Orthoplan Microscope - Included Items**

**Additional Items Included But Not Shown:**

- PN 10490 Hex Key 5/64 inch (for securing illuminator assembly to adapter)
- PN 10747 Hex Key 3mm long arm (for attaching adapter to microscope)



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		11/29/2021
PN 12117 Leitz Orthoplan Illuminator System		REV 3
SHEET 1		OF 1