


Nanodyne Replacement Illuminator for Nikon Alphaphot-2 Microscope: Included Items

Parts included but not shown:

- 8mm washer for PN 10213 pot
- 8mm nut for PN 10213 pot
- PN 10457 Hex Key 2.5mm short arm (for M3 socket head screws)
- PN 10337 Hex Key 2mm short arm (for M3 button head screws)



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Han-Seung Yang		11/20/2020
PN 11053 Nikon Alphaphot-2 Illuminator System		REV 4
SHEET 1		OF 1