


Nanodyne Replacement Illuminator for Nikon MM40 50W Microscope: Included Items

Additional Items Included But Not Shown:
 PN 10490 Hex Key 5/64 Short Arm (for set screw)



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.



PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

NOTE - You can specify Side Adjust when ordering to get the intensity adjust on the side of the illuminator, for easier installation.

Han-Seung Yang		10/8/2020
PN 11184 Nikon MM40 50W illuminator System - Remote pot	REV	3
SHEET 1 OF 1		