


Nanodyne Replacement Illuminator for Nikon Optiphot - Included Items



Additional Items Included But Not Shown:

- PN 10460 Hex Key 4mm (to remove bottom cover)
- PN 10457 Hex Key 2.5mm (for pot knob)
- PN 10490 Hex Key 5/64 (for 8-32 set screw)



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



3M VHB tape to secure pot cable

Note that this illuminator replaces the original 50W illuminator with the black metal housing. See our Nikon Optiphot-2 system if you have a 100W illuminator. That unit has a plastic housing that is usually cream colored, but some are black.

Han-Seung Yang		2/23/2021
PN 11223 Nikon Optiphot Illuminator System - Remote Pot		REV 5
SHEET 1		OF 1