


Nanodyne Replacement Illuminator for Nikon Epiphot: Included Items

Additional Items Included But Not Shown:

- PN 10399 Hex Key 5mm (to remove bottom cover)
- PN 10457 Hex Key 2.5mm (for pot knob)
- PN 10490 Hex Key 5/64 (for 8-32 set screw)



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Han-Seung Yang 2/22/2021	
PN 11478 Nikon Epiphot Illuminator System	REV 1
SHEET 1 OF 1	