


Nanodyne Replacement Illuminator for Nikon C-DSS Diascopic Microscope Base: Included Items

Additional Items Included But Not Shown:

- PN 10456 Hex Key 1.5mm (for M2 screws)
- PN 10457 Hex Key 2.5mm Short Arm (to remove OEM lamp housing)



The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Han-Seung Yang		9/3/2019
PN 11479 Nikon C-DSS Diascopic Illuminator System		REV 1
SHEET 1		OF 1