

Nanodyne Replacement Illuminator for Nikon E200 Microscope: Included Items

Additional Items Included But Not Shown:

- PN 10337 Hex Key 2mm Short Arm
- PN 10566 Hex Key 2.5mm long arm ball end

B


B

A

A

A

A



*** Extra tapes needed to secure pot plate.**

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



Tape to secure wires.

Han-Seung Yang 6/17/2020	
PN 11513 Nikon E200 Illuminator System	REV 1
SHEET 1 OF 1	