


# Nanodyne Replacement Illuminator for Nikon E800 Dual Illuminator Set Up - Included Items

## Additional Items Included But Not Shown:

- PN 10456 Hex Key 1.5mm (for pot knob)
- PN 10471 Hex Key 2mm x 9 inch T-Handle (for securing illuminator to scope)



Tape to secure wires.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10736 Rubber plug to block unused AC power receptacle.

Han-Seung Yang		11/16/2020
PN 11694 Nikon E800 Illuminator System - Dual Illuminator		REV 1
SHEET 1		OF 1