


# Nanodyne Replacement Illuminator for Nikon E800 - Included Items

## Additional Items Included But Not Shown:

- PN 10456 Hex Key 1.5mm (for pot knob)
- PN 10471 Hex Key 2mm x 9 inch T-Handle (for securing illuminator to scope)



Tape to secure wires.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10736 Rubber plug to block unused AC power receptacle.

PN 11243 Knob 0.236 in shaft 0.625 dia CLR Gloss

PN 11332 E800 Reflected Light Pot Plate

PN 11693 Nikon E800 illuminator Assy (HP3)

E800 illuminator assy has a diffuser while E600 illuminator assy doesn't.

PN 11170 Pot-Cable Assy - short shaft with tab (30 in cable)

Han-Seung Yang		3/31/2023	
PN 11694 Nikon E800 Illuminator System			REV 1
SHEET 1		OF 1	