


Nanodyne Replacement Illuminator for Nikon Epiphot 200 50W: Included Items

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 (for 8-32 set screw)



PN 11163 Power Supply 5V 2.1A and
PN 10734 Power Cable 1.35mm
ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Han-Seung Yang 2/22/2021	
PN 11783 Nikon Epiphot 200 50W Illuminator System	REV 1
SHEET 1 OF 1	