


**Nanodyne Replacement Illuminator for Nikon Ci-L Microscope: Included Items**

**Additional Items Included But Not Shown:**

- PN 10456 Hex Key 1.5mm (for pot knob and bottom cover)
- PN 10337 Hex Key 2mm (for bottom cover)
- PN 10457 Hex Key 2.5mm Short Arm (for M3 screws)



The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.



Tape to secure wires.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang 12/3/2020	
PN 11874 Nikon Ci-L Illuminator System	REV 1
SHEET 1 OF 1	