


Nanodyne Replacement Illuminator for Nikon Ti-U / Ti-S - Included Items

Additional Items Included But Not Shown:
 PN 10470 Hex Key 2mm Long Arm L



Tape to secure wires.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10736 Rubber plug to block unused AC power receptacle.

Han-Seung Yang		6/7/2021
PN 12057 Nikon Ti-U / Ti-S Illuminator System		REV 1
SHEET 1		OF 1