

**Nanodyne Replacement Illuminator for Nikon Diascopic Microscope illuminator System: Included Items**

**Additional Items Included But Not Shown:**

PN 10490 Hex Key 5/64 Short Arm (for set screw in illuminator - to adjust focus)

PN 10715 Olympus BH Illuminator Assembly (side adjustment)

PN 10467 Set Screw 8-32 x 0.25 brass tip

PN 11657 Set Screw Brass Tip 8-32 x 0.125

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

PN 12193 Nikon Diascopic Illuminator Adapter

PN 11656 Zeiss OPMI 1-FC Adapter Assy

Nikon Diascopic Stand OEM Port

Power supplies subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		4/20/2022
PN 12194 Nikon Diascopic Illuminator System		REV 1
SHEET 1		OF 1