


4 3 2 1

Nanodyne Replacement Illuminator for Nikon Optiphot Top - Included Items

PN 11145 T2PIHP3 Illuminator Assy
w Offset Lens - Side Adjust



PN 10743 SHCS M3 x 0.5 x 14mm (3 pcs)

PN 12070 Nikon Optiphot Adapter
for 100W Hg Lamp Bayonet

PN 10467 Set Screw
8-32 Brass Tip

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch (for securing illuminator to adapter)

PN 10566 Hex Key 2.5mm long arm ball end (for M3 SHCS to attach adapter to microscope)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and
PN 10734 Power Cable 1.35mm
ID/3.5mm OD x USB A 6 ft.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10736 Rubber plug to block
unused AC power receptacle.

Han-Seung Yang 3/28/2024	
PN 12536 Nikon Optiphot Top Illuminator System	REV 1
SHEET 1 OF 1	