


Nanodyne Replacement Illuminator for Olympus BX50-51-60-61 Top Illuminator: Included Items



Power supplies are subject to change due to availability and regulations.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.



PN 10736 Rubber plug to block unused AC power receptacle.

Wayne Bonin	5/23/2022
PN 10401 Olympus BX50-51-60-61 Top Illuminator System	REV 2
SHEET 1 OF 1	