

Nanodyne Replacement Illuminator for Olympus BX40 Microscope: Included Items



Tape to secure wires.




PN 10736 Rubber plug to block unused AC power receptacle.



PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.



PN 10908 BX40 Pot Plate - Cable Assy (9 inch long cable)

Wayne Bonin	3/11/2019
PN 10450 Olympus BX40 Illuminator System	REV 3
SHEET 1 OF 1	