

Nanodyne Replacement Illuminator for Olympus BX41/45 Microscope: Included Items

Additional Items Included But Not Shown:

PN 10334 Olympus T2PI
BX40 Series Illuminator Assy

PN 10456 Hex Key 1.5mm
3 pcs PN 10610 cable ties

***Note: Nanodyne's Illuminator System PN 10472 is used for both the Olympus BX41 and BX45.**

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

Pot Knob from BX41
(save from original microscope)

PN 11185 Custom Nut M7 x 0.75mm
- 9mm Bushing - 0.4375 in Hex

PN 10578 Serrated
Washer Assy



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.

PN 11170 Pot-Cable Assy
(shaft cut 5mm, 15 inch long)

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	6/11/2021
PN 10472 Olympus BX41 and 45 Illuminator system	REV 7