


Nanodyne Replacement Illuminator for Olympus BH-2 Microscope: Included Items



Additional Items Included But Not Shown:
 PN 10457 Hex Key 2.5mm (for M3 mounting screws)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.




PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

Power supplies subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Pot Knob from BH-2 (save from microscope)



PN 11002 Pot Cable Assy (15 in cable)

PN 10276 Slide Pot Assembly

PN 10274 M3 x 0.5 x 6mm (2 pcs) (Part of PN 10276)



PN 10736 Rubber plug to block unused AC power receptacle.

Wayne Bonin	3/12/2019
PN 10474 Olympus BH-2 Illuminator System	REV 5
SHEET 1	OF 1