


Nanodyne Replacement Illuminator for Olympus CHBS Microscope - Included Items

Additional Items Included But Not Shown:

PN 10457 Hex Key 2.5mm (for M3 mounting screws)



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Wayne Bonin	4/17/2019
PN 10491 Olympus CHBS Illuminator System	REV 5
SHEET 1 OF 1	