


Nanodyne Replacement Illuminator for Olympus CH30/CX40 Microscope: Included Items

Also Included but not shown:

PN 10457 Hex Key 2.5mm Short Arm (for M3 mounting screws)
 3M VHB Tape for securing cables after installation



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output



PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot



PN 10736 Rubber plug to block unused AC power receptacle.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	3/16/2019
PN 10556 Olympus CH30/CX40 Illuminator System	REV 3
SHEET 1 OF 1	