


# Nanodyne Replacement Illuminator for Olympus CH30/CX40 Microscope: Included Items

**Also Included but not shown:**  
 PN 10337 Hex Key 2mm Short Arm (for M3 BHCS)

**Tools not included but required:**  
 Phillips Screwdriver  
 11 mm socket  
 Misc tools for removing panels and original illuminator from microscope.



3M VHB Tape to secure wires.



The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

PN 10733 Power Supply & PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 ft.

Wayne Bonin	1/18/2020
PN 10556 Olympus CH30/CX40 Illuminator System	REV 4
SHEET 1 OF 1	