


Nanodyne Replacement Illuminator for Olympus BH Microscope - Included Items

Additional Items Included But Not Shown:
 PN 10457 Hex Key 2.5mm (for M3 mounting screws)



The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	11/24/2019
PN 10716 Olympus BH Illuminator System	REV 6
SHEET 1 OF 1	