

**Nanodyne Replacement Illuminator for Olympus BX43/46 Microscope (OEM adjust): Included Items**

**Additional Items Included But Not Shown:**

- PN 10490 Hex Key 5/64 inch (for set screw)
- PN 10566 Hex Key 2.5mm Long Arm (for M4 screw)
- PN 10747 Hex Key 3mm
- PN 10460 Hex Key 4mm

Note: diffuser in lens assy.

PN 11017 Olympus BX43 Illuminator Assy

PN 11015 Pot-Plate-Cable Assembly (15 in cable)

PN 12081 FHCS M4 x 0.7 x 10mm

Pot Knob (save from microscope)



Tape to secure wires.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	3/13/2023
PN 11018 Olympus BX43_46 Illuminator System (OEM Adjust)	REV 4
SHEET 1	OF 1