

Nanodyne Replacement Illuminator for Olympus BX43/46 Microscope (easy install): Included Items


Additional Items Included But Not Shown:

- PN 11721 Hex Key 1/8 inch (for set screw)
- PN 10747 Hex Key 3mm
- PN 10460 Hex Key 4mm


PN 11021 SHCS
M4 x 0.7 x 10mm





PN 11017
Olympus BX43_46
Illuminator Assy



PN 11054
Knob 0.50 dia.



PN 11062
Pot-Plate-Cable Assy
BX43_46 (easy install)
(15 in cable)



3M VHB Tape to secure wires.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	4/4/2020
PN 11018 Olympus BX43_46 Illuminator System (easy install)	REV 5
SHEET 1 OF 1	