

Nanodyne Replacement Illuminator for Olympus FH Microscope - Included Items



PN 11277 Olympus FH illuminator Assy



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	4/6/2019
PN 11275 Olympus FH Illuminator System	REV 2
SHEET 1	OF 1