


Nanodyne Replacement Illuminator for Olympus BH-2 Microscope Upper Illuminator (high power) - (15W OEM Lamp): Included Items



Additional Items Included But Not Shown:


- PN 10490 Hex Key 5/64 inch for set screw
- PN 10747 Hex Key 3mm long arm



PN 11163 Power Supply 5V 2.1A and
PN 10734 Power Cable 1.35mm
ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Han-Seung Yang	12/28/2020	
PN 11553 BH-2 Upper illuminator System (high power) - (15W OEM Lamp)		REV 1
SHEET 1 OF 1		