

Nanodyne Replacement Illuminator for Olympus CK40 Microscope - Included Items



PN 10736 Rubber plug to block unused AC power receptacle.




Tape to secure wires.



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.



PN 10334 Olympus BX40 Series Illuminator Assy

PN 11170 Pot-Cable Assy (40 inch long) - detabbed and full length shaft

Han-Seung Yang	3/12/2020
PN 11705 Olympus CK40 Illuminator System	REV 1
SHEET 1 OF 1	