

Nanodyne Replacement Illuminator for Olympus CKX41 Microscope - Included Items



PN 10736 Rubber plug to block unused AC power receptacle.



PCB light cover installed.




Tape to secure wires.




PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.




PN 10334 Olympus BX40 Series Illuminator Assy



PN 11170 Pot-Cable Assy (40 inch long) - detabbed and full length shaft

PN 12022 Olympus CKX41 PCB Light Cover-1

PN 12023 Olympus CKX41 PCB Light Cover-2



| | | |
|---|--|-----------|
| Han-Seung Yang | | 4/28/2021 |
| PN 11706 Olympus CKX41 Illuminator System | | REV 1 |
| SHEET 1 | | OF 1 |