


Nanodyne Replacement Illuminator for Olympus BX40 Microscope - Slide Pot: Included Items



PN 10334 Olympus BX40 Series Illuminator Assy



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.




PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.

PN 11002 Pot Cable Assy (12 in cable)



PN 11779 Slide Pot Assembly for BX40

Han-Seung Yang	6/26/2020
PN 11785 Olympus BX40 Illuminator System - Slide Pot	REV 1
SHEET 1 OF 1	