


**Nanodyne Replacement Illuminator for Olympus CHK-G Microscope - Included Items**

**Additional Items Included But Not Shown:**

PN 10457 Hex Key 2.5mm (for M3 mounting screws)



PN 12020 Olympus CHK-G Illuminator Assy

PN 10351 Pot-Cable Assy  
(shaft cut 5mm, detabbed)  
(8 inch long)



PN 10733 Power Supply - XP  
Power 5V 1A and PN 10734 Cable  
Assy 1.35mm ID x 3.5mm OD RA  
plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Han-Seung Yang 4/20/2021	
PN 12021 Olympus CHK-G Illuminator System	REV 1
SHEET 1 OF 1	