


Nanodyne Replacement Illuminator for Olympus BX53 Microscope - Halogen OEM Lamp (OEM adjust): Included Items

PN 10407 Olympus BX50-51 Bottom Illuminator Assy

Note: diffuser in lens assy.



Tape to secure wires.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	3/18/2022
PN 12186 Olympus BX53 Illuminator System - halogen OEM lamp (OEM Adjust)	REV 1
SHEET 1 OF 1	