


# Nanodyne Replacement Illuminator for Olympus IX70 Bottom Illuminator: Included Items

PN 10394 Olympus BX50-51 Top Illuminator Assy



3M VHB tape to secure pot cable

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID 3.5mm OD x USB A 6 ft.



PN 11795 Remote Control Box Assy (upon request)



PN 10736 Rubber plug to block unused AC power receptacle.

Power supplies are subject to change due to availability and regulations.

Han-Seung Yang		2/13/2023
PN 12195 Olympus IX70 Bottom Illuminator System		REV 1
SHEET 1		OF 1