


Nanodyne Replacement Illuminator for Olympus IX70 Top Microscope: Included Items



PN 10393 HP3 Illuminator Assy
Olympus BHS

Included Items Not Shown:

PN 10566 Hex Key 2.5mm long arm ball end (for 3 mm screws to attach OEM collector lens)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and
PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.



PN 11795 Remote Control Box Assy (42 inch)



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power receptacle.



PN 10275 SHCS M3 x 0.5 x 8mm (4 pcs)
(to attach OEM collector lens to illuminator)

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		2/9/2023
PN 12196 Olympus IX70 Top Illuminator System		REV 1
SHEET 1 OF 1		