

Nanodyne Replacement Illuminator for Olympus CKX41 Microscope Bottom High Power - Included Items

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch for set screw
PN 10566 Hex Key 2.5mm for M3 screw (to attach adapter to microscope)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 11163 Power Supply 5V 2.1A and
PN 10734 Power Cable 1.35mm
ID/3.5mm OD x USB A 6 ft.



Tape to secure wires.



PN 10736 Rubber plug to block
unused AC power receptacle.

Han-Seung Yang 3/16/2023	
PN 12379 Olympus CKX41 Bottom High Power Illuminator System	REV 1
SHEET 1 OF 1	

Power supplies are subject to substitution
without notice due to availability issues and
changes in regulations.