



Nanodyne Replacement Illuminator for StereoZoom Microscope Ring Light - Included Items



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output



PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot



Typical installation on a stereozoom microscope - microscope not included.

The illuminator may be powered by plugging the cable into the power supply shown above, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to supply issues and changes in regulations.

Wayne Bonin	12/22/2016
PN 10681 Ring Light Illuminator System	REV 1
SHEET 1 OF 1	