

3

 $\Delta$ 

2

## Additional Items Included But Not Shown:

- PN 10457 Hex Key 2.5mm (for M3 mounting screws)



В

K-

A

PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



Tape to secure wires.

Wayne Bonin	4/5/2019	9	
PN10494 Zeiss KF 2 Illuminator System			REV 4
SHEE	T 1 O	F	1